MULTI-FUNCTION INSTRUMENTS MODEL TA465 SERIES

The Airflow Model TA465 series are portable, handheld, Multi-Function Ventilation Test Instruments featuring a menu-driven user interface for easy operation in your local language. On-screen prompts and step-by-step instructions guide the user through instrument setup, operation, and field calibration. The TA465 also features an ergonomic, overmolded case design with probe holder and a keypad lockout to prevent tampering during unattended use. These instruments are available with or without a differential pressure sensor and are designed to work with a wide range of plug-in probes.




Features and Benefits

- + Best-in-class air velocity accuracy
- + Optional "smart" plug-in probes, including VOC, CO₂, and rotating vane probes
- + Accommodates up to two K-alloy thermocouples
- + Large graphic display
 - Displays up to five measurements simultaneously
 - On-screen messages and instructions
 - Program for local language
- + Intuitive menu structure allows for ease of use and setup
- + Multiple data logging formats
- + Bluetooth communications for transferring data or remote polling
- + Includes LogDat2 $^{\text{\tiny{M}}}$ downloading software with USB cable

Applications

- + HVAC testing and balancing
- + Clean room testing
- + Biological safety cabinet and laboratory fume hood testing
- + HVAC commissioning and troubleshooting
- + IAQ investigations
- + Thermal comfort studies
- + Ventilation evaluations
- + Process air flow testing





Airflow Model TA465 Plug-In Probes

The plug-in probes allow users to make various measurements by simply plugging in a different probe that has the features and functions best suited for a particular application.

Plug-in probes for the TA465 series can be ordered at any time and include a data sheet with certificate of traceability. When it's time for servicing, only the probe needs to be returned since all the calibration data is stored within the probe.

Thermoanemometer Air Velocity Probes

Airflow Instruments offers four models featuring multiple measurements in a compact, robust probe design. These telescopic probes are available in straight or articulating construction, and with or without a relative humidity sensor. Models with a relative humidity sensor can also calculate wet bulb and dewpoint temperature.

Common applications include duct traversing, face velocity testing of chemical fume hoods, biological safety cabinets and HEPA filters. When combined with the TA465, advanced measurement applications can be performed including heat flow, draft rate and turbulence intensity.

Rotating Vane Anemometer Probes

The 100 mm (4 in.) rotating vane probe measures air velocity and temperature with flow calculation. Measurement applications include face velocity as well as air velocity in turbulent airstreams. An optional telescopic articulating probe and an Aircone kit are also available.

Pitot Probes and Airflow Probe 800187

Pitot probes are used to obtain air velocity and air volume measurements within ductwork by performing a duct traverse. Consult factory for sizes and part numbers.

The Airflow Probe Model 800187 is an 46 cm (18 in.) straight Pitot probe that can be used to perform duct traverses and are ideally suited for measuring in small diameter ductwork.

LogDat2™ Downloading Software

The Airflow Model
TA465 Series includes
downloading software
called LogDat2. LogDat2
software transfers
the stored data from
the Model TA465
to a computer as a
spreadsheet file. This

Reading Type	Standard								
neading type	Temperature	70.0deg F							
	Pressure	29.92inHg							
Statistics	Channel:	Vel		Т		н		Dewpoint	Wetbulb
	Units:	ft/min		degF		%rh		degF	deg F
	Average:		827		71.9		22.1	31.3	51.7
	Minimum:		806		71.9		22.1	31.3	51.6
Date	Time	Vel		T		Н		Dewpoint	Wetbulb
MM/dd/yyyy	hh:mm:ss	ft/min		degF		%rh		deg F	deg F
3/1/2011	8:41:38	3	828		71.9		22.1	31.3	51.6
3/1/2011	8:41:40)	842		71.9		22.1	31.3	51.6
3/1/2011	8:41:42		836		71.9		22.1	31.3	51.6
3/1/2011	8:41:44		809		71.9		22.1	31.3	51.6
3/1/2011	8:41:46	5	806		71.9		22.1	31.3	51.6
3/1/2011	8:41:48	3	819		71.9		22.1	31.3	51.7
3/1/2011	8:41:50)	838		71.9		22.1	31.3	51.7
3/1/2011	8:41:52		837		71.9		22.2	31.3	51.3

software is useful for applications such as duct traverses, fume hood, and filter face velocity testing.

Data Collection and Reporting

Expanded data logging capacity and the inclusion of LogDat2 Downloading Software provides the capabilities to work more effectively and efficiently. The TA465 can store up to 38.9 days of data collected at one-minute log intervals. The stored data can be recalled, reviewed on screen, and downloaded for easy reporting.

- + Log multiple parameters to investigate trends.
- + Store up to 38.9 days of data collected at one-minute log intervals
- + User-selectable logging intervals
- + Download data to LogDat2 downloading software
- + Report generation



PROBE SPECIFICATIONS

960 Thermoanemometer Straight Probe Velocity and Temperature

Range 0 to 50 m/s (0 to 9,999 ft/min),

-18 to 93°C (0 to 200°F)

Accuracy ± 0.015 m/s ($\pm 3\%$ of reading or ± 3 ft/min),

whichever is greater485

±0.3°C (±0.5°F)6

Resolution 0.01 m/s (1 ft/min), 0.1°C (0.1°F)

962 Thermoanemometer Articulating Probe Velocity and Temperature

Range 0 to 50 m/s (0 to 9,999 ft/min),

-18 to 93°C (0 to 200°F)

Accuracy ± 0.015 m/s ($\pm 3\%$ of reading or ± 3 ft/min),

whichever is greater485

±0.3°C (±0.5°F)6

Resolution 0.01 m/s (1 ft/min), 0.1°C (0.1°F)

964 Thermoanemometer Straight Probe Velocity, Temperature and Humidity

Range 0 to 50 m/s (0 to 9,999 ft/min),

-10 to 60°C (14 to 140°F), 5 to 95% RH

Accuracy $\pm 3\%$ of reading or ± 0.015 m/s (± 3 ft/min),

whichever is greater485

±0.3°C (±0.5°F)⁶ +3% RH⁷

Resolution 0.01 m/s (1 ft/min), 0.1°C (0.1°F), 0.1% RH

966 Thermoanemometer Articulating Probe Velocity, Temperature and Humidity

Range 0 to 50 m/s (0 to 9,999 ft/min),

-10 to 60°C (14 to 140°F), 5 to 95% RH

Accuracy $\pm 3\%$ of reading or ± 0.015 m/s (± 3 ft/min),

whichever is greater,485

±0.3°C (±0.5°F),6 ±3% RH7

Resolution 0.01 m/s (1 ft/min), 0.1°F (0.1°C), 0.1%RH

995 Rotating Vane 4 in. (100 mm) Probe Velocity and Temperature $\,$

Range 0.25 to 30 m/s (50 to 6,000 ft/min),

0 to 60°C (32 to 140°F)

Accuracy $\pm 1\%$ of reading ± 0.02 m/s (± 4 ft/min),

±1.0°C (±2.0°F)

Resolution 0.01 m/s (1 ft/min), 0.1°C (0.1°F)

980 IAQ Probes CO₂, Temperature and Humidity

Range 0 to 5,000 ppm CO₂, 5 to 95% RH,

-10 to 60°C (14 to 140°F)

Accuracy $\pm 3\%$ of reading or ± 50 ppm CO₂,

whichever is greater9

±3% RH7

±0.5°C (±1.0°F)⁶

Resolution 1 ppm CO₂, 0.1°C (0.1°F), 0.1% RH

982 IAQ Probes Model CO, CO₂, Temperature and Humidity

Range $0 \text{ to } 500 \text{ ppm CO}, 0 \text{ to } 5000 \text{ ppm CO}_{2}$

-10 to 60°C (14 to 140°F), 5 to 95% RH

Accuracy $\pm 3\%$ of reading or ± 3 ppm CO,

whichever is greater8

±3% of reading or ±50 ppm CO₂,

whichever is greater⁹ ± 0.5 °C (± 1.0 °F)⁶ ± 3 % RH⁷

Resolution 0.1 ppm CO, 1 ppm CO₂, 0.1°C (0.1°F), 0.1% RH

792 and 794 Thermocouple Probes Temperature

Range -40 to 650°C (-40 to 1200°F) Accuracy ±0.056% of reading +1.1°C

 $(\pm 0.1\% \text{ of reading } + 2^{\circ}\text{F})$

Resolution 0.1°C (0.1°F)

984 Low Concentration (ppb) VOC and Temperature

Range 10 to 20,000 ppb, -10 to 60°C (14 to 140°F)

Accuracy ±0.5°C (±1.0°F)⁶ Resolution 10 ppb¹⁰, 0.1°C (0.1°F)

985 High Concentration (ppm) VOC and Temperature

Range 1 to 2,000 ppm, -10 to 60°C (14 to 140°F)

Accuracy ± 0.5 °C $(\pm 1.0$ °F)⁶ Resolution $1 \text{ ppm}^{10}, 0.1$ °C (0.1°F)

986 Low Concentration (ppb) VOC, Temperature, CO, and Humidity

Range 10 to 20,000 ppb VOC, 0 to 5,000 ppm CO₂,

-10 to 60°C (14 to 140°F), 5 to 95% RH

Accuracy $\pm 3\%$ of reading or 50 ppm $CO_{2'}$

whichever is greater ± 0.5 °C $(\pm 1.0$ °F)⁶ ± 3 % RH⁷

Resolution 10 ppb¹⁰ VOC, 0.1 ppm CO₂,

0.1°C (0.1°F), 0.1% RH

987 High Concentration (ppm) VOC, Temperature, CO, and Humidity

Range $1 \text{ to } 2,000 \text{ ppm VOC}, 0 \text{ to } 5,000 \text{ ppm CO}_{2'}$

-10 to 60°C (14 to 140°F), 5 to 95% RH

Accuracy $\pm 3\%$ of reading or 50 ppm CO₂,

whichever is greater $\pm 0.5^{\circ}\text{C} (\pm 1.0^{\circ}\text{F})^{6}$

Resolution 1 ppm 10 VOC, 0.1 ppm $^{CO}_2$,

0.1°C (0.1°F), 0.1% RH

SPECIFICATIONS

AIRFLOW MULTI-FUNCTION INSTRUMENTS MODELS TA465, TA465-A, TA465-P, TA465-X AND OPTIONAL PROBES

MODELS TA465, TA465-A, TA465-P, TA464-X AND OPTIONAL PROBES

Velocity (Pitot or Airflow probe for Meter Models 465, 465-A, 465-P)

Range¹ 1.27 to 78.7 m/s (250 to 15,500 ft/min) Accuracy² 10.16 m/s (±1.5% at 2,000 ft/min)

Resolution 0.01 m/s (1 ft/min)

Duct Size

Dimensions 2.5 to 1,270 cm in increments of 0.1 cm (1 to 500 inches in increments of 0.1 in.)

Volumetric Flow Rate

Actual range is a function of velocity, Range

pressure, duct size, and K factor

Static/Differential Pressure (Meter Models 465, 465-A, 465-P)

Range³ -28.0 to +28.0 mm Hg

-3,735 to +3,735 Pa (-15 to +15 in. H₂0)

Accuracy

 ± 0.01 mm Hg, ± 1 Pa ($\pm 1\%$ of reading ± 0.005 in. H₂0) Resolution 0.1 Pa, 0.01 mm Hg (0.001 in. H₂0)

Barometric Pressure

517.15 to 930.87 mm Hg Range (20.36 to 36.648 in. Hg)

Accuracy ±2% of reading

Instrument Temperature Range

Operating (Electronics) 5 to 45°C (40 to 113°F) Storage -20 to 60°C (-4 to 140°F)

Data Storage Capabilities

Range 26,500+ samples and 100 test IDs

Logging Interval

1 second to 1 hour

Time Constant

User selectable

External Meter Dimensions

9.7 cm x 21.1 cm x 5.3 cm (3.8 in. x 8.3 in. x 2.1 in.)

Meter Weight with Batteries

0.36 kg (0.8 lbs.)

Power Requirements

Four AA-size batteries or AC adapter

TO ORDER

Multi-Function Ventilation Meter with differential pressure sensor and Thermoanemometer Probe

Specify Description

TA465 Multi-function ventilation meter TA465-P with

straight air velocity probe Model 964

TA465-A Multi-function ventilation meter TA465-P with articulated air velocity probe Model 966

Multi-function Ventilation Meter Only. Choose a probe most appropriate for your measurement needs.

Specify Description

TA465-X Multi-function ventilation meter, no plug-in probes,

no differential pressure sensor

TA465-P Multi-function ventilation meter, no plug-in probes, with differential pressure sensor, tubing and static

NOTE: All models include: Instrument, hard carrying case, 4 alkaline batteries, USB cable, universal power supply, instruction manual, calibration certificate, LogDat2 downloading software.

Models TA465, TA465-A, and TA465-P also include (1) 2.4-m (8-ft) rubber tube and (1) static pressure tip.

- ¹ Pressure velocity measurements are not recommended below 5 m/s (1.000 ft/min) and are best suited to velocities over 10.00 m/s (2,000 ft/min). Range can vary depending on barometric pressure.
- ² Accuracy is a function of converting pressure to velocity. Conversion accuracy improves when actual pressure values increase.
- ³ Overpressure range = 360 mmHg, 48 kPa (190 in. H20).
- 4 Temperature compensated over an air temperature range of 5 to 65°C (40 to 150°F).
- 5 The accuracy statement begins at 0.15 m/s through 50 m/s
- (30 ft/min through 9,999 ft/min).
- ⁶ Accuracy with instrument case at 25°C (77°F), add uncertainty of 0.03°C/°C (0.05°F/°F) for change in instrument temperature.
- 7 Accuracy with probe at 25°C (77°F). Add uncertainty of 0.2% RH/°C (0.1% RH/°F) for change in probe temperature. Includes 1% hysteresis.
- ⁸ At 25°C (77°F). Add uncertainty of 0.36%/°C (±0.2%/°F) for change in temperature.
- 9 At calibration temperature, Add uncertainty of 0.5%/°C (±0.28%/°F) for change in temperature.
- 10 When response factor is set to 1.00.

Specifications subject to change without notice.

TSI and the TSI logo are registered trademarks, and Airflow, the Airflow logo and LogDat2 are trademarks of TSI Incorporated.



Airflow Instruments, TSI Instruments Ltd.

Visit our website at www.airflowinstruments.co.uk for more information.

Tel: +44 149 4 459200 Germany

Tel: +33 1 41 19 21 99

P/N 5001395 Rev D (A4)

©2015 TSI Incorporated

Tel: +49 241 523030